

Phenoxyacid Herbicides on Ultra Cyano (Overlay)

Peaks	RT (min)	Conc. (µg/mL)			
1. Dalapon	3.49	100	12. Dichlorprop	14.88	100
2. Picloram	5.98	100	13. 2,4,5-T	15.66	100
3. Dicamba	8.68	100	14. 2,4-DB	15.83	100
4. Chloramben	9.04	100	15. 2,4,5-TP	16.73	100
5. 4-Nitrophenol	9.71	100	16. Phenol, pentachloro	19.42	100
6. 2,4-Dichlorophenylacetic acid	12.53	40	17. Acifluorfen	19.47	100
7. Bentazon	12.70	100	18. Dinoseb	20.00	100
8. (4-Chloro-o-toloxyl)acetic acid	13.07	100	19. 2,4-D Butoxyethyl ester	20.78	200
9. 3,5-Dichlorobenzoic acid	13.41	100	20. 2,4,5-T Butoxyethyl ester	21.97	200
10. 2,4-D	13.49	100	21. 2,4,5-T Butyl ester	22.29	200
11. Mecoprop	14.68	100	22. 2,4-D 2-Ethylhexyl ester	24.12	200

